

BIOPHARMACEUTICAL FORMULATION



Spring 2020 | FTIR-PR03

ACCELERATE YOUR DEVELOPMENT OF STABLE PROTEIN FORMULATION WITH PROTA.

THE PROBLEM

Formulation of protein-based drugs often requires processes such as lyophilization or spray-drying. However, it has been shown that the freezing and drying steps involved in lyophilization could lead to protein unfolding. This unfolding of the native structure leads to an incomplete recovery of activity and results in aggregation and poor storage stability. These problems are often discerned through tedious and time-consuming measurements of activity or through detailed structure analysis. New and faster analytical techniques are needed in order to accelerate the development of protein formulation.

THE TECHNOLOGY

Infrared (IR) spectroscopy has widely proven to be a fast, cost-effective, accurate and reliable technique to determine secondary structure of a protein or to follow structural changes with perturbations like formulation. Amide I IR bands are a direct probe of the protein structure. Sample preparation is minimal and results can be obtained in less than half an hour with concentrations as low as 0.1-3mg/ml.

IR spectroscopy is the only spectroscopic technique that can be used for analysis of solids, such as lyophilized powders. Powder, fibrils, KBr disks, mulls or any other formulation can be used for study with FTIR. The use of the same technique provides a direct comparison between the native and the formulated protein.



THE SOLUTION

The PROTA Analyzer features dedicated software for protein analysis. Its data analysis function provides the necessary tools to quickly evaluate the structural differences between native and formulated proteins. This can be accomplished using the integrated formulation algorithm developed by Prof. John Carpenter, of the University of Colorado, or through the quantitative determination of structure using the built-in spectral database. This database includes the FT-IR spectra of 50 proteins with known secondary structure, and was acquired under strictly controlled conditions that match the crystallographic conditions. PROTA's database is the largest commercial protein IR database available.

BIOPHARMACEUTICAL FORMULATION



Spring 2020 | FTIR-PRO03

WE PROUDLY SERVE HUNDREDS OF COMPANIES, UNIVERSITIES AND GOVERNMENT AGENCIES

Over 5,000 determinations have been carried out over the last few years and many results have been submitted and accepted by legal and regulatory agencies. Various research has been published in top scientific publications.

AB BioTechnologies
Abbott
AbbVie
Akron BioTech
Alvogen. Inc
Amgen
Apex
AstraZeneca
Baxter
BioChem Pharma
Biogen-IDEC
Bristol-Myers Squibb
Cleveland Cord Blood Center
Complete Pharmaceuticals
Covance Biotechnology
Curagen
Duke University
DynPort
Eli Lilly
Endo Pharmaceuticals

Genzyme
GlaxoSmithKline
Halozyme
Hospira (A Pfizer Company)
Human Genome Sciences
Imclone
ImmunoGen, Inc.
InSight Biopharma
Ironwood Pharmaceuticals
KBI BioPharma
Labyrinth Biopharma
Magellan
Mars Inc.
MedImmune
Merck
Millenium Pharmaceuticals
Mission Barns
National Cancer Institute
Novartis
Novo Nordisk

Novotide
Oregon State University
Pfizer
Plastic Surgery Innovations
Regeneron
Royal Danish School of Pharmacy
Silk Technologies
Shire Biologics
Teva
Trimeris
University of Alabama
University of Colorado
University of Kansas
University of N. Colorado
University of Texas Medical Center
USDA-NCAUR
Wyeth (now Pfizer)
Xeris Pharmaceuticals
Zymogenetics

We specialize in the analysis of...

proteins, peptides, carbohydrates, DNA, and chiral molecules by VCD, protein FT-IR, Optical Rotation, and UV Visible Spectroscopy.

BioTools offers **lab services** for absolute configuration of chiral molecules. You provide the sample and structure, we do the measurement and calculation. Receive a detailed report for submission to regulatory agencies, patent protection, publication or internal results.

Visit www.biotoools.us to get started!

IT'S AS EASY AS...

1

Fill out an analysis form w/ structure

2

Mail in your samples

3

Receive a detailed signed report



Characterization Experts: Chirality & Biologics

PROUDLY SERVING THE SCIENTIFIC COMMUNITY FOR OVER 20 YEARS!

17546 Bee Line Highway | Jupiter, Florida 33458 USA
ph. 561.625.0133 | fax. 561.625.0717 | info@biotoools.us | www.biotoools.us