


ChiralIR-2XTM

Vibrational Circular Dichroism Spectrometer


- #1 Selling VCD instrument *In The World*
- #1 Choice of Pharma & Academics
- #1 Published VCD Instrument
- Thousands of molecules measured using the ChiralIRTM platform



Discover VCD with the ChiralIR-2XTM and look at molecules with a new *twist* of light!

Unique Advantages of the ChiralIR-2X™

- Dual Source increases signal-to-noise ratio (S/N) of approximately 3.6 times over any single-source system
- Exclusive stage for elimination of cell artifacts with SyncRoCell™
- Complete spectral range of measurement 4000-850 cm^{-1} (no change of filters required)
- 30-48 hour hold detector time
- Simultaneous collection of IR & VCD spectra (AC & DC)



Sample spectrum in full range